

Real-world Measures of Functional Gait in Patients with Amyotrophic Lateral Sclerosis Using Wearable Digital Health Technologies



¹Rakesh Pilkar, PhD, ²Cory Holdom, ¹Dawid Gerstel, PhD, ¹Matthew R Patterson, PhD, ¹Ali Neishabouri, PhD¹, ¹Andy Liu, PhD, ¹Sylvain Zorman, PhD, ²Frederik Steyn, PhD, ¹Christine Guo, PhD,

¹ActiGraph L.LC, Pensacola, FL 32502, USA,

²The University of Queensland, Brisbane, Queensland, Australia

Introduction

- Amyotrophic lateral sclerosis (ALS) is a neurodegenerative disease that progressively damages parts of the nervous system, causes muscle degradation, and significantly affects movement, speech and respiratory functions.
- Quantification of gait, an essential component for understanding ALS progression as well as intervention outcomes, is usually performed by subjective clinical evaluations or burdensome laboratory-based experiments.
- Wearable digital health technologies (DHTs) can reduce the patient burden of lab visits and provide objective and sensitive continuous data to quantify real-world behavior.





Fig. 1. The Voice of the Patient, ALS drug dev survey, the ALS association, Oct, 2019. **A.** Patients' response, **B**. Caregivers' response

To demonstrate the feasibility of unobtrusive wrist-based Objective accelerometry to free-living gait metrics in persons with ALS (pwALS).

- Q1: Can wrist-based digital gait measures distinguish between ALS and the controls?
- Q2: Can digital gait measures track the ALS progression?
- Q3: How well the digital gait measures correlate with the clinical measure, ALS-FRS?

Methods

Setting: Free-living ambulation

Device: GT9X (ActiGraph, L.L.C., Pensacola, FL)

Procedure:


- GT9X activity monitors were placed on non-dominant wrists of 104 pwALS and 54 agematched controls.
- The accelerometer data were continuously recorded for 7 to 8 days at 30 Hz.

Group	n	age (years)	height (cm)	weight (kg)
control	54	55.62 (11.7)	170.57 (8.17)	79.45 (20.96
ALS	104	60.14 (8.53)	173.9 (8.3)	79.11 (17.75

Data processing and outcome measures


Accelerometer data were processed using a custom-built algorithm to extract step count, gait speed, and cadence.

Fig. 2. Data processing pipelines for deriving PA and gait measures.




Results


Q1: Can wrist-based digital gait measures distinguish between ALS and the controls?




Q3: How to the digital gait measures correlate with the clinical measure?



Q2: Can digital gait measures track the ALS progression?





Conclusions & future work

- Wrist-accelerometry can provide rich data to extract objective and sensitive measures of real-world gait for pwALS
 - ✓ Showed significant differences between the control and pwALS
 - ✓ Can track ALS progression
 - ✓ Significantly correlated with the clinical measure ALS-FRS
 - Bout-definition can influence the variability and effect size for gait metrics (speed, cadence) in the freeliving condition. To be investigated.



- 1. The Voice of the Patient, ALS drug dev survey, the ALS association, Oct, 2019
- 2. Holdom CJ, van Unnik JWJ, van Eijk RPA, van den Berg LH, Henderson RD, Ngo ST, Steyn FJ. Use of hip- versus wrist-based actigraphy for assessing functional decline and disease progression in patients with motor neuron disease. J Neurol. 2023 Feb 6. doi: 10.1007/s00415-023-11584-7. Epub ahead 3. R. Delgado-Gonzalo et al., "Physical activity profiling: Activity-specific step counting and energy expenditure models using 3D wrist
- acceleration," 2015 37th Annual International Conference of the IEEE Engineering in Medicine and Biology Society (EMBC), Milan, Italy, 2015, pp.
- 4. M. Bertschi et al., "Accurate walking and running speed estimation using wrist inertial data," 2015 37th Annual International Conference of the IEEE Engineering in Medicine and Biology Society (EMBC), Milan, Italy, 2015, pp. 8083-8086, doi: 10.1109/EMBC.2015.7320269.